

coOrdination®


The system of coOrdination® Titanium double sleeves indicates titanium sleeves for axial guiding of cylindrical drills in surgical templates or surgical guides. Under certain conditions they might be used together with depth controlled surgical instruments within 3D implant planning tools.

Therefore the exact drill diameter and the length between drill tip to depth stop have to be known. The distance between implant shoulder and the drill sleeve is specified by the free drill length, the sleeve length and the implant length.

To use double sleeves with depth stop drills the diameter of the depth stop has to be at least 5.0 mm. Otherwise it would stop inside the sleeves funnel. It has to be made sure that the drill fits the sleeve correctly in advance. Conical shaped drills cannot be guided in a cylindrical sleeve.

Requirements

- o Which drill shall be used? Does the drill fit the sleeve? (check in advance!)
- o Which sleeve should be used?
- o Does the drill have a depth stop > 5.0 mm?
- o Distance between drill tip and depth stop is larger than implant plus sleeve length.



By using the template drill for inner or outer sleeve a hole is made the gives the shape of the sleeves outer surface. The sleeves upper collar has a height of 0.2 mm. By choosing the sleeve diameter within the software it has to be considered, that the inner sleeve sits 0.45 mm on the outer sleeve. Depending on the usage of inner sleeve and/or outer sleeve this distance has to be considerate within planning. All provided measures may vary due to production tolerances.

Axis adjustment

The gonyX™ axis adjustment is independent of vertical sleeve position and can be used from the software directly. The screenshot to the right shows an example for the outer sleeve.



<p>coDiagnostiX™ Patientdata</p> <p>Version: 8.0 Name: Steco Licensed to: Date of birth: - Institut Straumann AG (150... Patient ID: 123456789</p>			
<p>Plan: Oberkiefer Position: 36</p> <p>Sleeve: steco system-technik Titan Outer-Sleeves Artikelnummer: M.27.02.D350 Sleeve length: 6.0mm Diameter (inner): 4.0mm (3.50mm)</p> <p>Depthcontrol: Zero level to implant tip: 30.7mm Drilling depth for sleeve: 11.7mm Check value (top of sleeve): 5.7mm Sleeve to bone (space foringiva): 7.0mm Zero level to implant tip: 30.7mm</p> <p>Implant: Straumann SLActive Bone Level (Regular CrossFit™) Artikelnummer: 021.6112/G Length: 12.0mm Diameter 1: 4.8mm Diameter 2: 4.8mm</p> <p>Surgical protocol: Drilling depth to stop: 25.0mm</p> <p>Reference coordinates: X: 18.4mm; Y: 18.6mm; Alpha: 76.7°; Beta: 7.7°</p>		<p>gonyX™ templateplan FD Inotation (World Dental Federation)</p> <p>D: 76.7° 90° 60°</p> <p>C: 7.7° 20° 10° 0°</p> <p>B: 205.4° 210° 180°</p> <p>A: 14.5° 30° 0°</p> <p>Done: <input type="checkbox"/> Verified: <input type="checkbox"/> Editor: _____</p> <p style="font-size: small;">Copyright © 2011 Institut Straumann AG. All rights reserved. Printed: 2012-02-27 19:38</p>	

Titanium double sleeve with depth stop in Straumann® coDiagnostiX™

The following statement apply to the common use of outer and inner sleeves.

- The gonyX™ axis adjustment is independent of vertical sleeve position and can be used from the software directly.
- Take bearing of zero level with the tip of template drill (M.27.02.B400)!
- Drill depth is calculated like following:

check value upper sleeve rim from software
 + 7.2 mm* (tip of template drill to end of bevel)
 + 0.2 mm collar of outer sleeve
 + 0.45 mm collar of inner sleeve**
 = drill depth



* This length definition is only sketchy because the measure is taken within the spiral cut of the drill. The real measure is depending on the drill procedure as well. Make sure to plan additional safety distance!



** Inner sleeve bears on outer sleeve.

Titanium inner sleeve with depth stop in Straumann® coDiagnostiX™

The following statement apply to the use of inner sleeves without outer sleeves.

- The gonyX™ axis adjustment is independent of vertical sleeve position and can be used from the software directly.
- Take bearing of zero level with the tip of template drill (M.27.03.B350)!
- Drill depth is calculated like following:

check value upper sleeve rim from software
 + 10.8 mm* (tip of template drill to end of bevel)
 + 0.2 mm collar of inner sleeve**
 = drill depth



* This length definition is only sketchy because the measure is taken within the spiral cut of the drill. The real measure is depending on the drill procedure as well. Make sure to plan additional safety distance!

** Inner sleeve bears on drilled hole.